

Rascal Rueba

Mid-range powerchairs



Your new route to independence

RASCAL 

Rascal Rueba

Rueba Specification:

Type Class	B / 2
Overall length	107cm / 42.1"
Overall width	61cm / 24"
Floor to seat height	43 - 45.5cm (min-max)
Seat angle setting	0 - 15 °
Seat width	45.5cm / 18"
Seat depth	56cm / 22"
Maximum carrying capacity	159kgs / 25 stone
Battery voltage & capacity	12V / 55Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	7cm / 2.8"
Maximum kerb height	5cm / 2"
Range*	up to 33.5km / 20.8 miles
The maximum safe slope	6°
Turn-around width	1.60m / 5'2"
Colour	Silver Blue

*Range on full charge and flat ground

Rascal Rueba Profile:

The Rascal Rueba 4 mph Rear-Wheel Drive powerchair suitable for indoor and outdoor use comes with 55Ah batteries as standard and features a flexible seating system with a removable 18" seat cushion enabling bespoke adaptations. The capacity to adjust the seat tilt, backrest angle and the horizontal strap tension for total comfort achieves the ideal seating position. The narrow and stream-lined design fits easily through doorways offering full independence.

Local Rascal Mobility Retailer:

Main Features:

- Adjustable upholstery on back
- Removable armrests
- Coil spring suspension
- Convenient joystick charging port
- Pneumatic tyres & anti-tip castors
- Back angle adjustable from 0 - 20 ° in 5 ° increments
- Manual adjustment seat tilt 2 - 18 °
- Seat folds for easy transportation
- Swing-away legrests

**OPTIONAL
ACCESSORY**

Optional Extras:

- Elevating legrests
- Attendant control bracket
- Lighting kit
- Capability to fit 16"-22" seat cushion



Specifications subject to change without prior notice
Images are for illustrative purposes only
E. and O.E- v3 Final-July 2023

electricmobility.co.uk