

# VelumiLi



Lithium powered scooters



**LiFetime Guarantee**

on Lithium Ferro-Phosphate Battery\*




The Lithium Ferro-Phosphate Battery is **'built to last'** & the **SAFEST** Lithium Battery in the Mobility sector

**RASCAL**



# Rascal VelumiLi Profile:

The **Rascal VelumiLi**  is a lightweight, Aluminium boot scooter with high quality State-Of-The-Art **Lithium Ferro-Phosphate Battery**. At just 28.8kg\*, our lightest Rascal ever also includes a LiFetime Battery Guarantee for peace of mind. Brimming with features the slimline frame is just 18.6" wide yet feels sturdy with 8" dual compound tyres providing a smooth ride, comfortable & spacious seating position for travelling & stability for everyday adventures. Discover your style with the choice of 6 vibrant colours including the **New Silver Ice** option.

Featuring two easy access charging points, on the tiller or on the removable battery box, which is ideal for indoor charging. **The compact & clever styling with super light front section weighing just 10.2kg** makes it exceptional to manoeuvre, with 1150mm turning radius perfect for travelling near or far, at home or overseas.


## Enjoy the convenience & impressive durability...




**...Today is the day to start your journey!**

Introducing  
**Silver Ice**  
12.4Ah



With the **Rascal VelumiLi**  every detail makes a statement, from the crafted white stitching on the seat, LED front light to be seen with, reflective side strip & a good size front basket which is quick to attach on or off for bitesize essentials. This 4mph lightweight scooter has a range of upto 13 miles\*\*, ideal for both indoor & outdoor use, it is slick & easy to dismantle, with a T Bar tiller locking feature for efficient transportation. It may be small but it is mighty & great value for money.

The  **BRITISH MADE** Lithium Ferro-Phosphate Battery is **'built to last'** & the **SAFEST** Lithium Battery in the Mobility sector

\* Weight excluding Battery † LiFetime Battery Guarantee T&Cs apply \*\* 12.4Ah-13 miles / 10Ah-11.9 miles / 25Ah-26.2 miles

## Main Features:

- Lightweight, aluminium, high quality rugged frame
- British Made / Medical /Military grade **Lithium Ferro-Phosphate Battery**
- Airline Compliant 12.4Ah battery weighs just 4.1kgs / 9 lbs
- 1150mm turning radius for exceptional manoeuvrability
- Tiller locking for easy transportation
- Fold down backrest & simple to remove adjustable armrests
- Front basket has good size & depth
- Two easy access charging points on tiller & battery box
- Dual compound puncture proof tyres for impressive stability
- 15.3" seat height that offers real comfort & support for long journeys
- Folded Dimensions (LxWxH) 106cm x 47cm x 71cm / 41.7" x 18.6" x 27.9"



Clever, practical **front handle** increases accessibility



## Air travel & CAA Lithium Battery Legislation

Airline Compliant 10Ah & 12.4Ah battery options meet I.A.T.A. UN38.3 Passenger Aircraft Transportation requirements for overseas travel & **can be taken on a plane.** (Subject to Airline Regulation & Restrictions)

## 10Ah features:

- Airline Compliant 10Ah battery weighs just 3.7kgs / 7.7 lbs
- Range<sup>^</sup> upto 19.2km / 11.9 miles
- Product weight including battery 32.3kgs / 71.2 lbs

## 25Ah features:

- 25Ah NOT Airline Compliant
- Range<sup>^</sup> upto 42.1km / 26.2 miles
- Product weight including battery 36.7kgs / 80.9 lbs

<sup>^</sup> Range on full charge and flat ground



Safely transport your belongings & travel essentials with this compact accessory

OPTIONAL




Extras:

- Rascal arm rest bag



## Boot Transportable

The **Rascal VelumiLi**  with high quality State-Of-The-Art **Lithium Ferro-Phosphate** (LiFePO<sub>4</sub>) Battery. Relax & enjoy the ride in the knowledge that the **Lithium Ferro-Phosphate** (LiFePO<sub>4</sub>) Battery pack is the safest, **greenest** with an increased efficiency & performance. Plus the LiFetime Battery Guarantee†.

## Your perfect travel companion



Silver Ice 10Ah



Specification:	10Ah	12.4Ah	25Ah
Type Class	B / 2		
Overall length	106cm / 41.7"		
Overall width	47cm / 18.6"		
Overall height	88cm / 34.7"		
Floor to seat height	39 - 43cm (min-max)		
Seat width	41cm / 16.1"		
Seat depth	37cm / 14.6"		
Folded Dimensions LxWxH	106cm / 47cm / 71cm (41.7" / 18.6" / 27.9")		
Maximum carrying capacity	115kgs / 18 stone		
Battery voltage and capacity	25.6V / 10Ah	25.6V / 6.2Ah x 2	25.6V / 25Ah
Maximum speed	6.4km/h / 4mph		
Ground clearance	4.5cm / 1.8"		
Maximum kerb height	5cm / 2"		
Range*	up to 19.2km / 11.9miles	up to 21km / 13miles	up to 42.1km / 26.2miles
The maximum safe slope	6°		
Turn-around width	1.25m / 4'1"		
Colour	Blue Moon / Green Lightning / Orange Sunset / Pink Thunder / Yellow Star		
Single Heaviest Component (Front Section)	10.2kgs / 22.5 lbs		
Product weight (Exc battery) (Inc battery)	28.8kgs / 63.5 lbs 32.3kgs / 71.2 lbs	32.8kgs / 72.3 lbs	36.7kgs / 80.9 lbs

\* Range on full charge and flat ground

See full Tech Page on Website & Owner's Manual for details



LiFetime Battery Guarantee

## Today is the day to begin your adventure

Local Rascal Retailer:



† LiFetime Battery Guarantee T&Cs apply