


Rascal Veo X



Rascal Veo X

Rascal Veo X Specification:

Type Class	B / 2
Overall length	108cm / 42.5"
Overall width	60cm / 23.6"
Floor to seat height	47-49cm (min-max)
Seat width	45.7cm / 18"
Seat depth	40.6cm / 16"
Maximum carrying capacity	129.2kgs / 20 stone
Battery voltage and capacity	12V / 22Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	6cm / 2.3"
Maximum kerb height	5cm / 2"
Range*	up to 21km / 13.5 miles
The maximum safe slope	8°
Turn-around width	1.40m / 4'6"
Colour	Red / Teal

*Range on full charge and flat ground.

Main Features:

- Car transportable
- Fold down backrest
- Flip-up width adjustable armrests
- Removable front basket
- Tiller adjustment
- Solid tyres all-round
- Delta-bar tiller
- Single heaviest component (front chassis) 17.4kgs / 38.3lbs



Rascal Veo X Profile:

The Rascal Veo X is a robust and highly manoeuvrable travel scooter that is great for both indoor and outdoor use. It has among the best legroom of any scooter in its class and a really easy to use 'take apart' mechanism. With its larger battery capacity giving up to a 13.5 mile range, upgraded seat and wheels, delta-bar tiller and striking modern styling it is ideal for users who need more range and a little more comfort from their ordinary car transportable scooter. The Veo X now features new wide dual density wheels for an improved ride.



Local Rascal Mobility Retailer:

E. and O.E- v 5 Final-July 2015