

Rascal 388 XL



Mid-size scooters

Your new route to independence



Rascal 388 XL

Rascal 388 XL Specification:

Type Class	C / 3
Overall length	125cm / 49"
Overall width	54cm / 21"
Seat width	45.7cm / 18"
Seat depth	45.7cm / 18"
Maximum carrying capacity	136kgs / 21 stone
Battery voltage and capacity	12V / 33Ah x 2
Maximum speed	9.65km/h / 6mph
Ground clearance	13cm / 5"
Maximum kerb height	5cm / 2"
Range*	up to 32km / 20 miles
The maximum safe slope	8°
Turn-around width	2.87m / 9'5"
Colour	Black / Red

*Range on full charge and flat ground.

Main Features:


- Outstanding reliability
- Flip-up width adjustable armrests
- Comfortable high back, height adjustable, sliding swivel seat
- Narrow compact chassis
- Lift-up tiller head for easy access
- Upgraded captains seat
- Smooth controls
- One rear view mirror as standard
- Front basket as standard



Rascal 388 XL Profile:

The Rascal 388 XL is the perfect Class 3 scooter with 6 mph top speed, stylish 11" wheels, full road lighting and with its narrow chassis it combines safety with performance and manoeuvrability. Sturdier in build and ideal for people requiring additional comfort on those slightly longer journeys - a superb scooter for days out and leisure trips. With a comfortable high back seat with headrest, width adjustable armrests and tiller, this scooter can be 'adapted' to suit the individual's requirements. 388 S and 388 Deluxe models also available.

Local Rascal Mobility Retailer:



Specifications subject to change without prior notice
E. and O.E - v 5 Final-November 2017