

# Rascal Vantage X



**Mid-size scooters**

**Your new route to independence**

# Rascal Vantage X

## Rascal Vantage X Specification:

Type Class	C / 3
Overall length	120cm / 47.2"
Overall width	58.5cm / 23"
Seat width	45.7cm / 18"
Seat depth	40.6cm / 16"
Maximum carrying capacity	149.7kgs / 23.5 stone
Battery voltage and capacity	12V / 50Ah x 2
Maximum speed	9.65km/h / 6mph
Ground clearance	6.5cm / 2.6"
Maximum kerb height	7.6cm / 3"
Range*	up to 42km / 26 miles
The maximum safe slope	10°
Turn-around width	1.40m / 4'6"
Colour	Black / Red

\*Range on full charge and flat ground.

## Main Features:

- Fold down backrest
- Flip-up width adjustable armrests
- Removable front basket
- Infinite adjustable delta tiller position
- Modern delta handlebars
- Waterproof panel
- Emergency hydraulic braking system

## Optional Extras:

- Rain canopy
- Cane and crutch holder bag



## Rascal Vantage X Profile:

The Rascal Vantage X is a compact, value for money 6 mph road scooter with extremely smooth user-friendly controls. With a comfortable captain's seat as standard, infinitely adjustable delta-bar tiller steering and waterproof switches this highly manoeuvrable scooter features LED lighting and can adequately accommodate users up to 23.5 stone.



Local Rascal Mobility Retailer:

Specifications subject to change without prior notice  
E. and O.E.- v 5 Final-October 2017