

Rascal Frontier



Performance scooters

Your new route to independence

Rascal Frontier

Rascal Frontier Specification:

Type Class	C / 3
Overall length	143cm / 56"
Overall width	65cm / 26"
Seat width	45.7cm / 18"
Seat depth	45.7cm / 18"
Maximum carrying capacity	136kgs / 21 stone
Battery voltage and capacity**	12V / 33Ah x 2
Maximum speed	12.9km/h / 8mph
Ground clearance	9cm / 3.5"
Maximum kerb height	10cm / 4"
Range*- standard battery pack (33Ah)	up to 24km / 15 miles
Range*- upgraded battery pack (55Ah)	up to 40km / 25 miles
The maximum safe slope	10°
Turn-around width	1.70m / 5'6"
Colour	Black / Red

*Range on full charge and flat ground. **Upgrade to 12V / 55Ah x 2

Main Features:

- Flip-up width adjustable armrests
- Leaf spring suspension
- Removable front basket
- Articulated tiller adjustment with useful battery charging port
- Modern delta handlebars
- Waterproof panel and twin mirrors
- Large shock absorbing tyres for maximum comfort on the move

Optional Extras:

- A stylish rear mounted lockable box fitted on a sliding mount to enable easy access without adjusting the seat
- Standard rain canopy - upgrade to a polycarbonate front shield for an additional cost



Rascal Frontier Profile:

The Rascal Frontier is an entry level 8 mph scooter, a proven performer with excellent legroom, comfort and dependability. It features leaf spring suspension, a delta-bar tiller with a comfortable seat as well as waterproof switches and a large basket making the Frontier one of the most popular scooters on the market. Available with an optional battery upgrade to 55Ah, that can give up to a 25 mile range.

Local Rascal Mobility Retailer:

Specifications subject to change without prior notice
E. and O.E- v 6 Final October 2017

