

# Ultralite 765



**Compact Powerchairs**

**Your new route to independence**

# Ultralite 765

## Rascal Ultralite 765 Specification:

Type Class	B / 2
Overall length	96cm / 38"
Overall width	57cm / 22"
Floor to seat height	50cm (min)
Seat width	45.7cm / 18"
Seat depth	43.1cm / 17"
Maximum carrying capacity	91kgs / 14 stone
Battery voltage and capacity	12V / 12Ah x 2
Maximum speed	6km/h / 3.8mph
Ground clearance	4cm / 1.5"
Range*	up to 14.4km / 9 miles
The maximum safe slope	6°
Colour	Red

\*Range on full charge and flat ground.


## Main Features:

- Car transportable
- Adjustable backrest
- Flip-up width adjustable armrests
- Handy storage pouch behind the seat
- Convenient joystick charging port
- Easy to charge batteries off-board
- Adjustable captains seat
- Flip-up footplate



## Ultralite 765 Profile:

The Ultralite 765 is a compact and lightweight powerchair that gives you the lifestyle you desire. Featuring easy to use controls, a simple 2-step folding process for storage and transportation and Rear-Wheel Drive giving superb manoeuvrability. The 765 powerchair gives great comfort with its high back seating, total reliability and easy positioning under tables and desks.



Local Rascal Mobility Retailer:

Specifications subject to change without prior notice  
E. and O.E- v 2 Final-August 2014