

Rascal P321



Compact powerchairs

Your new route to independence

Rascal P321

Rascal P321 Specification:

Type Class	B / 2
Overall length	88cm / 35"
Overall width	58.5cm / 23"
Seat width	45.7cm / 18"
Seat depth	40.6cm / 16"
Maximum carrying capacity	113kgs / 18 stone
Battery voltage and capacity	12V / 14Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	2.3cm / 0.9"
Maximum kerb height	2.5cm / 1"
Range*	up to 13km / 8 miles
The maximum safe slope	6°
Turn-around width	1.06m / 3'5"
Colour	Red

*Range on full charge and flat ground.


Main Features:

- Car transportable
- Flip-up fully adjustable armrests
- Solid tyres all round
- Convenient joystick charging port
- Folding seat
- Flip-up footplate



Rascal P321 Profile:

The Rascal P321 gives you a lightweight, compact Rear-Wheel Drive 4 mph powerchair. Specifically designed to combine both easy dismantling and portable mobility, the P321 is the ideal choice for indoor and occasional outdoor use and delivers versatility and total comfort. Getting the seat position right for you is simple with armrests, seat height and joystick adjustments.



Local Rascal Mobility Retailer:

Specifications subject to change without prior notice
E. and O.E - v 3 Final October 2017