

Rascal P327 Mini



Modern
Colour Accents




Up to 35km

22 Mile
Range*



Max. Carrying
Capacity
136kgs / 21 Stone



Full Suspension

Rascal P327 Mini

Rascal P327 Mini Specification:

Type Class	B / 2
Overall length	109cm / 42.9"
Overall width	62.4cm / 24.5"
Floor to seat height	52cm (min)
Seat width	45.7cm / 18"
Seat depth	43.1cm / 17"
Maximum carrying capacity	136kgs / 21 stone
Battery voltage and capacity	12V / 34Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	5cm / 2"
Maximum kerb height	5cm / 2"
Range*	up to 35km / 22 miles
The maximum safe slope	6°
Turn-around width	1.18m / 3'8"
Colour	Red / Teal

*Range on full charge and flat ground.

Main Features:

- Flip-up fully adjustable armrests
- Full suspension
- Convenient joystick charging port
- Dual density, puncture resistant drive wheels for maximum comfort on the move
- Swing-away joystick
- High back adjustable seat



Rascal P327 Mini Profile:

The Rascal P327 Mini is a compact 4 mph Mid-Wheel Drive indoor/outdoor powerchair, with full suspension and comfortable scooter style seating. The P327 Mini is the ideal choice for indoor and outdoor use where space is restricted. It delivers great performance with its 34Ah batteries contributing to its outstanding manoeuvrability and a smooth ride.



Local Rascal Mobility Retailer:

Specifications subject to change without prior notice
E. and O.E- v 1 Final-August 2014