


Also Available

Matching
Sofa

Also Available

Static
Armchair

Kensey

MADE TO
ORDER



Shown in
Waterfall Back

Made
in the UK





CosiChair

Kensey Size Guide:

| | Small | Medium | Large |
|----------------|-------------------|---------------------|---------------------|
| Seat height | 45.7cm / 18" | 50.8cm / 20" | 53.3cm / 21" |
| Seat width | 45.7cm / 18" | 50.8cm / 20" | 55.9cm / 22" |
| Seat depth | 45.7cm / 18" | 50.8cm / 20" | 53.3cm / 21" |
| Arm height | 20.3cm / 8" | 19.1cm / 7½" | 19.1cm / 7½" |
| Back height | 66cm / 26" | 68.6cm / 27" | 71.1cm / 28" |
| Overall width | 76.2cm / 30" | 81.3cm / 32" | 86.4cm / 34" |
| Overall height | 108cm / 42½" | 110.5cm / 43½" | 113cm / 44½" |
| Weight limit | 20 stone / 127kgs | 25 stone / 158.8kgs | 25 stone / 158.8kgs |

The Kensey is our new UK Premium Tilt-in-Space (TIS) Single-motor Cosi Chair providing a smooth rise and recline movement. It's Waterfall Back design (shown on the front) delivers support to the user's back, shoulders and neck areas throughout travel. Choose from any other back design below:



Main Features:

- Premium Seminar TIS mechanism with Variable Angle Lift
- Arm covers and Antimacassar included Free of Charge as standard
- Choice of any back design or style of fabric on a Made to order basis **only**
- 5-button power on handset
- New pleated wing design
- Heat & massage / drop in seat cushions for pressure care issues - both available on request*
- 5 Year Guarantee on frames, mechanisms and all electrics

Local Cosi Chair Retailer:

*Available on Made to order **only** - please ask for more details.
Please note colours may vary slightly from the images shown.