

Also Available



Matching Sofa

Also Available



Static Armchair

# Jubilee

Lateral & Waterfall Back



Shown in  
Lateral Back

Made  
in the UK



# Jubilee

## Lateral & Waterfall Back

Choose our furniture with confidence -  
All riser recliner chairs are tested by FIRA -  
The leading authority in furniture safety by  
design.



### Jubilee Size Guide:


\* Also available in Petite - measurements Seat H 16" x Seat W 18" x Seat D 16" - on Made to order

	Small	Medium	Large
Seat height	45.7cm / 18"	50.8cm / 20"	53.3cm / 21"
Seat width	45.7cm / 18"	50.8cm / 20"	55.9cm / 22"
Seat depth	45.7cm / 18"	50.8cm / 20"	53.3cm / 21"
Arm height	19.1cm / 7½"	19.1cm / 7½"	19.1cm / 7½"
Back height	68.6cm / 27"	73.7cm / 29"	78.7cm / 31"
Overall width	76.2cm / 30"	81.3cm / 32"	86.4cm / 34"
Overall height	108cm / 42½"	110.5cm / 43½"	113cm / 44½"
Weight limit	20 stone / 127kgs	25 stone / 158.8kgs	25 stone / 158.8kgs

The Jubilee is our new UK Premium Tilt-in-Space (TIS) Dual-motor Cosi Chair designed to maximise support throughout motion with a range of stylish back designs & upholstery options to provide exceptional postural seating positions to suit all individuals.

### Main Features:

- Premium Seminar TIS mechanism with Variable Angle Lift
- Arm covers and Antimacassar included Free of Charge as standard
- Choice of any back design or style of fabric on a Made to order basis **only**\*
- 5-button power on handset
- New pleated wing design
- Heat & massage / drop in seat cushions for pressure care issues - both available on request†
- 5 Year Guarantee on frames, mechanisms and all electrics



\*Button & Mini-Lateral Backs available on Made to order **only** - please ask for more details.



Local Cosi Chair Retailer:

†Available on Made to order **only** - please ask for more details.  
Please note colours may vary slightly from the images shown.



Electric Mobility Euro Ltd, Canal Way, Ilminster, Somerset TA19 9DL

0800 252614 sales@electricmobility.co.uk electricmobility.co.uk

Due to a policy of continual improvements, Electric Mobility Euro Ltd reserve the right to change specification without prior notice.

