

Banwell



Banwell

Choose our furniture with confidence -
All riser recliner chairs are tested by FIRA -
The leading authority in furniture safety by
design.



The Banwell is a powerful bariatric riser recliner chair for the heavier or wider customer that provides the upmost reassurance & comfort.

40 stone Banwell with the 25 stone leg lift – option 1

This has been designed with the larger customer in mind, those suffering from heavier legs that are swollen or Oedema and may require a footboard suitable to support the extra carrying capacity. Unlike other standard riser recliner chairs available, which can only carry approximately 20% total lifting capacity on the footboard, the CosiChair Banwell carries up to 25 stone on its footboard area. It is a dual motor chair controlling the base and back with a further 2 motors powering the leg lift only.

35 stone Banwell with Tilt-in-space – option 2

This dual motor recliner is for the heavier customer who requires unparalleled comfort & support of the head and back area. As you recline the chair, the seat to back angle doesn't alter allowing for full support throughout the chairs motion, it then allows the back to be separately adjusted for the ultimate comfort and positioning.



*Available to order in Seasons Cocoa. Other fabric options available on a Made to order **basis** - please ask for more details.

Banwell Size Guide for stocked product:

	40 stone option 1	35 stone option 2
Seat height	48.2cm / 19"	48.2cm / 19"
Seat width	63.5cm / 25"	63.5cm / 25"
Seat depth	53.3cm / 21"	53.3cm / 21"
Arm height	20.3cm / 8"	20.3cm / 8"
Back height	88.9cm / 35"	88.9cm / 35"
Overall width	81.3cm / 32"	81.3cm / 32"
Overall height	122cm / 48"	122cm / 48"
Weight limit	40 st / 254kgs	35 st / 222kgs

Main Features:

- 2 variations of mechanism available:
35 stone Tilt-in-Space or 40 stone with 25 stone leg lift
- Arm covers and Antimacassar included Free of Charge as standard
- Waterfall Back design with wings
- Made to order & Made to Measure available
- High Density Foam seat base
- Pressure care cushions available on request*
- 5 Year Guarantee on frames, mechanisms and all electrics

Local Cosi Chair Retailer:



Electric Mobility Euro Ltd, Canal Way, Ilminster, Somerset TA19 9DL

0800 252614 sales@electricmobility.co.uk electricmobility.co.uk

Due to a policy of continual improvements, Electric Mobility Euro Ltd reserve the right to change specification without prior notice.



MADE TO ORDER

