

# Rascal P327 SL



Max. Carrying  
Capacity  
182kgs / 28 Stone

**RASCAL**  
MOBILITY



**Mid-wheel powerchairs**

**Your new route to independence**

# Rascal P327 SL

## Rascal P327 SL Specification:

Type Class	B / 2
Overall length	112cm / 44'
Overall width	69cm / 27"
Floor to seat height	
Lowest to highest elevation position	62-87cm (min-max)
Seat width	50.8cm / 20"
Seat depth	45.7cm / 18"
Maximum carrying capacity	182kgs / 28 stone
Battery voltage and capacity	12V / 55Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	5cm / 2"
Maximum kerb height	5cm / 2'
Range*	up to 32km / 20 miles
The maximum safe slope	6°
Turn-around width	1.20 / 4'
Colour	Red / Teal

\*Range on full charge and flat ground.

## Main Features:

- Flip-up fully adjustable armrests
- Full suspension
- Swing-away joystick
- Convenient joystick charging port
- Adjustable captains seat
- Flip-up footplate
- **Maximum height when fully elevated - 147cm / 57.9"**



## Rascal P327 SL Profile:

The Rascal P327 SL is a heavy-duty 4 mph Mid-Wheel Drive powerchair suitable for indoor/outdoor usage, with durable full suspension. The SL has a range of up to 20mph and delivers great performance with 55Ah batteries. This model includes a comfortable captain seat with scissor lift mechanism to allow for exceptional manoeuvrability. The scissor lift delivers smooth and quiet lift action and reaches a maximum height of 147 cm when fully elevated.



Local Rascal Mobility Retailer:

Specifications subject to change without prior notice  
E. and O.E- v 1 Final - October 2017