

Rascal Ventura



Performance scooters

Your new route to independence

Rascal Ventura


Rascal Ventura Specification:

Type Class	C / 3
Overall length	131cm / 51.6"
Overall width	68.4cm / 26.9"
Seat width	45.7cm / 18"
Seat depth	40.6cm / 16"
Max carrying capacity - batt pk (50Ah)	159kgs / 25 stone
Battery voltage and capacity** †	12V / 50Ah x 2
Maximum speed	12.9km/h / 8mph
Ground clearance	9cm / 3.5"
Maximum kerb height	10cm / 4"
Range* - standard battery pack (50Ah)	up to 37km / 23 miles
Range* - upgraded battery pack (75/80Ah)	up to 59km / 36 miles
The maximum safe slope	10°
Turn-around width	1.85m / 6'
Colour	Black / Silver



*Range on full charge and flat ground. **Upgrade to 12V / 75/80Ah x 2

†Max carrying capacity - batt pk (75/80Ah) 133.8kgs / 21 stone

Main Features:

- Front coil and rear mono-shock suspension for a superior ride
- Gas strut tiller with easy access battery charging port
- Modern delta handlebars
- Modern LCD panel and twin mirrors
- Safety turning sensor, slows down scooter when turning at full speed
- Auxiliary USB power outlet for charging mobile phones etc 
- Emergency braking system
- Modern LED lighting system

Optional Extras:

- A stylish rear mounted lockable box fitted on a sliding mount to enable easy access without adjusting the seat 
- Standard rain canopy - upgrade to a polycarbonate front shield for an additional cost 



Rascal Ventura Profile:

The Rascal Ventura is a robust, comfortable 8 mph scooter that will take you wherever you wish to go. With its modern practical design and state of the art styling the Ventura is a real head turner. Couple this with its delta handlebar, tiller adjustment, LED lighting, LCD instrument panel and a waterproof USB power outlet the Ventura has all the makings of a contemporary classic.



Local Rascal Mobility Retailer:

Specifications subject to change without prior notice
E. and O.E. - v 6 Final-October 2017

