

Rascal P327 Seat Lift




Max. Carrying
Capacity
182kgs / 28 Stone

RASCAL
MOBILITY

Mid-wheel powerchairs

Your new route to independence

Rascal P327 Seat Lift

Rascal P327 Seat Lift Specification:

Type Class	B / 2
Overall length	112cm / 44'
Overall width	69cm / 27"
Floor to seat height	62-87cm (min-max)
Lowest to highest elevation position	
Seat width	51cm / 20"
Seat depth	46cm / 18"
Maximum carrying capacity	182kgs / 28 stone
Battery voltage and capacity	12V / 55Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	5cm / 2"
Maximum kerb height	5cm / 2"
Range*	up to 32km / 20 miles
The maximum safe slope	6°
Turn-around width	1.20m / 4'
Colour	Red / Teal

*Range on full charge and flat ground.

Main Features:

- Flip-up fully adjustable armrests
- Full suspension
- Swing-away joystick
- Convenient joystick charging port
- Adjustable Captain's seat
- Flip-up footplate
- Maximum floor to seat height when fully elevated - 87 cm / 34"



Rascal P327 Seat Lift (SL) Profile:

The Rascal P327 SL is a heavy-duty 4 mph Mid-Wheel Drive powerchair suitable for indoor/outdoor usage, with durable full suspension & outstanding manoeuvrability. The SL has a range of up to 20 miles and delivers great performance with 55Ah batteries. This model includes a comfortable Captain's seat with scissor lift mechanism. The scissor lift delivers smooth and quiet lift action and reaches a maximum floor to seat height of 87 cm when fully elevated.



Local Rascal Mobility Retailer:

Specifications subject to change without prior notice
E. and O.E- v 2 Final - December 2017