

# Breydon

CosiChair 



Made  
in the UK



# Breydon

Our static chairs are certified by FIRA - the leading authority on furniture testing. You can purchase with confidence knowing that our chairs have successfully passed the relevant strength, stability and durability tests.



## Breydon Size Guide:

	One size
Seat height	48.2cm / 19"
Seat width	49cm / 19¼"
Seat depth	48.2cm / 19"
Arm height	18.4cm / 7¼"
Back height	73.6cm / 29"
Overall width	59.7cm / 23½"
Overall height	119.3cm / 47"
Weight limit	18 stone / 114.3kgs



The Breydon fireside chair provides a low cost multi-function high seat chair for those who need the comfort of a robust high seat with carefully shaped arm cushions and long handgrips to assist in standing. The removable seat cushion enables the seating to be adapted to suit individuals who require specific pressure care solutions.

## Main Features:

- Solid seat base with a removable zipped cushion
- Exposed wooden arms for easy grip
- Shaped back for lower lumbar support with wings
- Solid beech frame with mortise and tenon joints offering strength and stability
- Other fabrics are available



Herringbone  
Cocoa



Herringbone  
Oyster



Solid Beech  
Frame



Foot Stool  
available

Please note colours may vary slightly from the images shown.

Local Cosi Chair Retailer:



Electric Mobility Ltd, Canal Way, Ilminster, Somerset, TA19 9DL

☎ 0800 252614 ✉ sales@electricmobility.co.uk ⓘ www.electricmobility.co.uk

📘 www.facebook.com/CosiChairs 🐦 www.twitter.com/CosiChairs

Due to a policy of continual improvements, Electric Mobility reserve the right to change specification without prior notice.

