

Also
Available



Matching
Sofa

Ambassador

Waterfall Back



Ambassador

Waterfall Back

The Ambassador offers outstanding levels of comfort and styling with its superbly flexible seating positioning system and deep generous waterfall cushions providing the upmost essential riser recliner chair. Its unique design reaches ultimate positioning by achieving the Trendelenburg position.

Ambassador Waterfall Back Size Guide:

	Small	Medium	Large
Seat height	45.7cm / 18"	50.8cm / 20"	53.3cm / 21"
Seat width	48.2cm / 19"	53.3cm / 21"	58.4cm / 23"
Seat depth	43.2cm / 17"	50.8cm / 20"	53.3cm / 21"
Arm height	16.5cm / 6½"	16.5cm / 6½"	15.2cm / 6"
Back height	66cm / 26"	68.6cm / 27"	71.1cm / 28"
Overall width	77.5cm / 30½"	81.9cm / 32¼"	86.4cm / 34"
Overall height	110.5cm / 43½"	114.3cm / 45"	119.3cm / 47"
Weight limit	21 stone / 136kgs	26.8 stone / 170kgs	26.8 stone / 170kgs

Main Features:

- Multi-functional Dual-motor 6-button handset, allows independent controls of the backrest and footrest
- Button or Waterfall Back design
- Canadian kiln-dried maple hardwood frame
- Heat and massage available on request
- One Year Guarantee on frames, mechanisms and all electrics including handsets
- Arm covers, Antimacassar and 2 pockets come as standard
- Trendelenburg position which gives high leg elevation and raises the feet higher than the heart
- 2-Seater static sofa available on Made to order



Local Cosi Chair Retailer:

Please note colours may vary slightly from the images shown.
The Ambassador is available from STOCK in all the above fabrics.



Electric Mobility Euro Ltd, Canal Way, Ilminster, Somerset TA19 9DL

0800 252614 sales@electricmobility.co.uk electricmobility.co.uk

Due to a policy of continual improvements, Electric Mobility Euro Ltd reserve the right to change specification without prior notice.

