

Also
Available



Matching
Sofa

Ellen



Made
in the UK




CosiChair

Ellen Size Guide:


	One size
Seat height	48.3cm / 19"
Seat width	50.8cm / 20"
Seat depth	50.8cm / 20"
Arm height	20.3cm / 8"
Back height	71.1cm / 28"
Overall width	82.5cm / 32½"
Overall height	111.8cm / 44"
Weight limit	20 stone / 127kgs



The Ellen Cosi Chair is a comfortable Tilt-in-Space Single-motor chair with pocket sprung chaise seating and a deep fibre filled Waterfall backrest. The design retains the same angle between the seat and backrest to maximise comfort for individuals with delicate skin conditions or where legs need elevating above the hips.

Main Features:

- 3-button power on handset
- Waterfall Back design
- Pocket sprung seat base
- Scroll arms with wings
- 5 Year Guarantee on frames, mechanisms and all electrics including handsets
- Arm covers and Antimacassar as standard
- 2-Seater static sofa available on Made to order
- Other fabrics are available on Made to order




Spray
Cocoa



Spray
Evergreen



Spray
Jute



Spray
Rose



Easy to use
handset



Recline
position

Please note colours may vary slightly from the images shown.

Local Cosi Chair Retailer: