## Rascal Rialto

Crash Test Approved for Forward Facing Applications: (ISO 7176-19) as an occupant seat in an appropriate Road Vehicle



Your versatile Powerchair

RASCAL

# Rascal Rialto

### **Rialto Specification:**



Type Class	B/2
Overall length	109cm / 42.9"
Overall width	62.4cm / 24.5"
Floor to seat height	52cm / 20.4" (min)
Seat width	46cm / 18.1" & 51cm / 20.1"
Seat depth	39.5cm / 15.6", 42cm / 16.5" & 45.5cm / 17.9"
Maximum carrying capacity	136kgs / 21.4 stone
Battery voltage and capacity	12V / 35Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	5cm / 2"
Maximum kerb height	5cm / 2"
Range*	up to 26.3km / 16.3 miles
The maximum safe slope	6° - Do not exceed

Lightning
Single Heaviest component 62kg / 136.7lbs

\*Range on full charge and flat ground

Turning circle

Colour

See full Tech Page on Website & Owner's Manual for details

#### Rascal Rialto Profile:

The Rascal Rialto is our proven, manoeuvrable indoors & outdoors 4mph Mid-Wheel Drive Powerchair offering configurable options normally found on higher cost models. The Rialto is an excellent personal transport proposition bridging the gap between high-end & basic powerchairs. The Rialto boasts Crash Test Approval & multi-size seat, width, depth & angle adjustment seating options.

1.18m / 3'87"

Sunset Orange/ Green





• Crash Tested: (ISO 7176-19)

Koller™ Restraint

- Configurable seating system
- Includes 18.75" seat cushion
- Swing-away detachable legrests
- Tension adjustable back
- Coil-sprung suspension
- Puncture proof tyres & castors
- Swing-away joystick
- Adjustable seat depth, back & tilt angle

#### OPTIONAL

- · Centre footplate optio
- · Attendant control

### OPTIONAL ACCESSORY

Armrest bag





Crash Test Approved\* for Forward Facing Applications: (ISO 7176-19)

Allows the use of the powerchair as a seat in a motor-vehicle when used in conjuction with an approved restraint system.

\* Technical T's & C's apply





