

# Liteway 8 Max Grip™



**Transportable scooters**

**Your new route to independence**

# Liteway 8 Max Grip™


## Liteway 8 Max Grip™ Specification:

Type Class	C / 3
Overall length	126cm / 49.6"
Overall width	60cm / 23.6"
Floor to seat height	63cm (min)
Seat width	45.7cm / 18"
Seat depth	45.7cm / 18"
Maximum carrying capacity	136kgs / 21 stone
Battery voltage and capacity	3 sets of 12V 14Ah x 2
Maximum speed	12.8km/h / 8mph
Ground clearance	8.5cm / 3.3"
Maximum kerb height	5cm / 2"
Range*	up to 38km / 23.8 miles
The maximum safe slope	10°
Turn-around width	1.98m / 6'5"
Colour	Blue / Graphite

\*Range on full charge and flat ground.

## Max Grip™ LSD - Limited Slip Differential Gearbox Features:

- Improved grip, handling & safety
- Unique take-apart construction & fully road legal
- Improved energy efficiency - increased range



## Main Features:

- Car Transportable / Road Legal
- Rear coil spring suspension
- High efficiency LED lighting
- Lightweight aluminium frame
- Fold-down tiller for easy storage & transportation
- Luxury seating with swivel function & flip up / width adjustable armrests
- Single heaviest component (rear chassis) 18.1kgs / 40lbs

## Liteway 8 Max Grip™ Profile:

At last, a Class 3, transportable, lightweight scooter that has been especially designed to go the extra mile. The Liteway 8 Max Grip™ provides versatility of a transportable scooter with the performance & range of a full Class 3 (road legal) unit normally too heavy & cumbersome to transport.

## Optional Extras:

- Front basket
- Super compact seat
- Off-board charging station



OPTIONAL  
ACCESSORY  
OR  
UPGRADE



Local Rascal Mobility Retailer:

Specifications subject to change without prior notice  
E. and O.E. v 4 Final-August 2017