

# Rascal Veo



Transportable scooters



Your new route to independence



# Rascal Veo

## Rascal Veo Specification:

Type Class	B / 2
Overall length	108cm / 42.5"
Overall width	50.5cm / 19.9"
Floor to seat height	47-49cm (min-max)
Seat width	43.1cm / 17"
Seat depth	35.5cm / 14"
Maximum carrying capacity	113kgs / 18 stone
Battery voltage and capacity	12V / 14Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	6cm / 2.3"
Maximum kerb height	5cm / 2"
Range*	up to 18.5km / 11.5 miles
The maximum safe slope	8°
Turn-around width	1.40m / 4'6"
Colour	Citrus / Orange / Red / Teal

\*Range on full charge and flat ground.

## Main Features:

- Car transportable
- Fold down backrest
- Flip-up width adjustable armrests
- Removable front basket
- Tiller adjustment
- Solid tyres all-round
- Delta-bar tiller
- Single heaviest component (front chassis) 17.4kgs / 38.3lbs



## Rascal Veo Profile:

The Rascal Veo is a robust and highly manoeuvrable travel scooter that is great for both indoor and outdoor use. It has among the best legroom of any scooter in its class and an easy to use 'take apart' mechanism. With its delta-bar tiller and striking modern styling, the Veo is the ideal solution if you regularly need to transport your scooter into the boot of a car.

Local Rascal Mobility Retailer:

Specifications subject to change without prior notice  
E. and O.E - v 6 Final-April 2018