

Rascal WeGo 250



Compact powerchairs

Your new route to independence

Rascal WeGo 250

Rascal WeGo 250 Specification:

Type Class	B / 2
Overall length	100cm / 39.4"
Overall width	60cm / 23.6"
Floor to seat height	50-58cm (min-max)
Seat width	45.7cm / 18"
Seat depth	40.6cm / 16"
Maximum carrying capacity	113kgs / 18 stone
Battery voltage and capacity	12V / 14Ah x 2
Maximum speed	4.8km/h / 3mph
Ground clearance	3.5cm / 1.4"
Maximum kerb height	2.5cm / 1"
Range*	up to 8.8km / 5.5 miles
The maximum safe slope	6°
Turn-around width	1.32 / 4'3"
Colour	Red

*Range on full charge and flat ground.

Main Features:

- Car transportable
- Flip-up width adjustable armrests
- Under-seat storage facility
- Easy to charge batteries off-board
- Solid tyres all round
- Flip-up footplate



Optional Extras:

- Optional high back Captain's seat
- Optional larger footplate



Rascal WeGo 250 Profile:

Pushing a manual wheelchair can be tiring especially on hills. Powered by its battery pack the Rascal WeGo 250 puts you in complete control as the delta-style handlebars allow your carer to manoeuvre with ease compared to a traditional style attendant electric wheelchair. Offering total security with its additional safety brake and being easy to dismantle and transport in the boot of a car, the WeGo 250 sits uniquely on its own in today's market.



Local Rascal Mobility Retailer:

Specifications subject to change without prior notice
E. and O.E- v 4 Final-September 2018

