

Rascal WeGo 250



Compact powerchairs

Your new route to independence

Rascal WeGo 250

Rascal WeGo 250 Specification:

| | |
|------------------------------|-------------------------|
| Type Class | B / 2 |
| Overall length | 100cm / 39.4" |
| Overall width | 60cm / 23.6" |
| Floor to seat height | 50-58cm (min-max) |
| Seat width | 45.7cm / 18" |
| Seat depth | 40.6cm / 16" |
| Maximum carrying capacity | 114kgs / 18 stone |
| Battery voltage and capacity | 12V / 15Ah x 2 |
| Maximum speed | 4.8km/h / 3mph |
| Ground clearance | 3.5cm / 1.4" |
| Maximum kerb height | 2.5cm / 1" |
| Range* | up to 8.8km / 5.5 miles |
| The maximum safe slope | 6° |
| Turn-around width | 1.32 / 4'3" |
| Colour | Red |

*Range on full charge and flat ground.

Main Features:

- Car transportable
- Flip-up width adjustable armrests
- Under-seat storage facility
- Easy to charge batteries off-board
- Solid tyres all round
- Flip-up footplate


Optional Extras:

- Optional high back captains seat
- Optional larger footplate



Rascal WeGo 250 Profile:

Pushing a manual wheelchair can be tiring especially on hills. Powered by its battery pack the Rascal WeGo 250 puts you in complete control as the delta-style handlebars allow the carer to manoeuvre with ease compared to a traditional style attendant electric wheelchair. Offering total security with its additional safety brake and being easy to dismantle and transport in the boot of a car, the WeGo 250 sits uniquely on its own in today's market.



Local Rascal Mobility Retailer:

Specifications subject to change without prior notice
E. and O.E- v 3 Final-May 2018

