

Rascal P312 Turnabout with Seat Lift



RASCAL
MOBILITY

Mid-range powerchairs

Your new route to independence

Rascal P312 Turnabout with Seat Lift

Rascal P312 Turnabout with Seat Lift Specification:

Type Class	B / 2
Overall length	96cm / 38"
Overall width	61cm / 24"
Floor to seat height with seat lift	57- 69.5cm (min-max)
Seat width	45.7cm / 18"
Seat depth	40.6cm / 16"
Maximum carrying capacity	159kgs / 25 stone
Battery voltage and capacity	12V / 50Ah x 2
Maximum speed	6km/h / 3.8mph
Ground clearance	5cm / 2"
Maximum kerb height	5cm / 2"
Range*	up to 29km / 18 miles
The maximum safe slope	10°
Turn-around width	1.18m / 3'10"
Colour	Red

*Range on full charge and flat ground.

Main Features:

- 5 inch Column Seat Lift
- Front-Wheel or Rear-Wheel Drive
- Flip-up width adjustable armrests
- Convenient joystick charging port
- Swing-away joystick
- High back adjustable seat
- Fold-away footplates front and rear

Optional Extras:

- Optional accessory swing-away legrests
- Optional accessory elevating legrests



Rascal P312 Turnabout with Seat Lift Profile:

The Rascal P312 Turnabout with Seat Lift gives flexibility to the user including a dual directional drive with a rotating seat giving you the option of Front-Wheel Drive for superb manoeuvrability indoors or Rear-Wheel Drive for extra traction outdoors. The P312 Turnabout Seat Lift offers well defined driving positions with a swing-away joystick and a fully adjustable seat and headrest for day-long comfort.



Local Rascal Mobility Retailer:

Specifications subject to change without prior notice
E. and O.E - v 4b Final-September 2017