

Rascal P327 XL



Full Suspension



Up to 32km


20 Mile
Range*



Max. Carrying
Capacity
204kgs / 32 Stone



Great
Turn-around
Width
1.20m / 4'



Rascal P327 XL

Rascal P327 XL Specification:

Type Class	B / 2
Overall length	112cm / 44"
Overall width	69cm / 27"
Floor to seat height	54-56.5cm (min-max)
Seat width	52cm / 20.5"
Seat depth	45.7cm / 18"
Maximum carrying capacity	204kgs / 32 stone
Battery voltage and capacity	12V / 55Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	5cm / 2"
Maximum kerb height	5cm / 2"
Range*	up to 32km / 20 miles
The maximum safe slope	6°
Turn-around width	1.20 / 4'
Colour	Red / Teal

*Range on full charge and flat ground.

Main Features:

- Flip-up fully adjustable armrests
- Full suspension
- Swing-away joystick
- Convenient joystick charging port
- Adjustable Captain's seat
- Flip-up footplate

Optional Extras:

- Optional accessory swing-away legrests
- Optional accessory elevating legrests
- Larger seat available - 56cm / 22"



Rascal P327 XL Profile:

The Rascal P327 XL is a heavy-duty 4 mph Mid-Wheel Drive outdoor/indoor powerchair, with durable full suspension and comfortable cushioned Captain's seat mounted on a larger base. The P327 XL delivers great performance with 55Ah batteries contributing to its excellent range, outstanding manoeuvrability and a smooth ride. P327 Mini, P327 (standard) and Seat Lift models also available.



Local Rascal Mobility Retailer:

Specifications subject to change without prior notice
E. and O.E- v 5 Final-September 2018