

# Rascal P200



**Mid-range powerchairs**

**Your new route to independence**

# Rascal P200

## Rascal P200 Specification:

Type Class	B / 2
Overall length	109cm / 42.9"
Overall width	65.1cm / 25.6" - 73.7cm / 29" (min - max)
Floor to seat height	57cm (min)
Seat width	45.7cm / 18"
Seat depth	40.6cm / 16"
Maximum carrying capacity	136kgs / 21 stone
Battery voltage and capacity	12V / 55Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	5cm / 2"
Maximum kerb height	5cm / 2"
Range*	up to 26km / 16.2 miles
The maximum safe slope	8°
Turn-around width	1.32m / 4'3"
Colour	Silver

\*Range on full charge and flat ground.

## Main Features:

- Adjustable upholstery on back
- Removable armrests
- Coil spring suspension
- Convenient joystick charging port
- Solid tyres all round
- Manual angle adjustable seat tilt
- Swing-away legrests

## Optional Extras:

- Optional accessory elevating legrests
- Option of an Attendant control bracket
- Option for Kerb Climbers




## Rascal P200 Profile:

The Rascal P200 4 mph Rear-Wheel Drive powerchair suitable for indoor and outdoor use comes with 55Ah batteries as standard and features a flexible seating system with a removable seat cushion enabling bespoke adaptations. The capacity to adjust the seat tilt, backrest angle and the horizontal strap tension for total comfort achieves the ideal seating position.

### Coil Suspension




Up to 26km



16.2 Mile  
Range\*

Local Rascal Mobility Retailer:

### Great Turn-around Width



### Swing-away Legrests



Specifications subject to change without prior notice  
E. and O.E- v 3b Final-September 2018

OPTIONAL  
ACCESSORY