

Rascal Veo Sport



RASCAL Exclusive
All-round suspension



RASCAL
MOBILITY

Transportable scooters

Your new route to independence

Rascal Veo Sport

Rascal Veo Sport Specification:

Type Class	B / 2
Overall length	108cm / 42.5"
Overall width	60cm / 23.6"
Floor to seat height	47-49cm (min-max)
Seat width	45.7cm / 18"
Seat depth	40.6cm / 16"
Maximum carrying capacity	129.2kgs / 20 stone
Battery voltage and capacity	12V / 22Ah or 14Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	6cm / 2.3"
Maximum kerb height	5cm / 2"
Range* - lightweight battery pack (14Ah)	up to 18.5km / 11.5 miles
Range* - standard battery pack (22Ah)	up to 21km / 13.5 miles
The maximum safe slope	8°
Turn-around width	1.40m / 4'6"
Colour	Aquamarine / Platinum Grey

*Range on full charge and flat ground.

Main Features:

- A RASCAL  Exclusive
- Car transportable
- All-round suspension
- Range up to 13.5 miles
- Stylish sports seat
- Dual compound puncture proof tyres
- Comfortable delta tiller
- Rear safety reflectors
- Colours in latest high-tech, two-tone matt finishes
- Optional upgrade - front lighting system
- **Now available with lightweight 14Ah battery pack**



Rascal Veo Sport Profile:

It's all about the drive.

The Veo Sport is fully loaded with features. A modern, sporty design with class leading all-round suspension. With 22Ah batteries, a superb capacity for a scooter of its size, affording a range of up to 13.5 miles. A padded delta tiller and luxury seat, detailed with unique white accent stitching, the Veo Sport is a real pleasure to ride. The Veo Sport's easy and robust take apart mechanism allows you to store or drive away with it in the boot of your car, so it is very practical yet with its distinctive looks, exclusive all-round suspension, **this is the travel scooter with the edge!**



Local Rascal Mobility Retailer:

Specifications subject to change without prior notice
E. and O.E- v 2 Final-August 2017