

Also Available



Matching Sofa

Also Available



Static Armchair

Kensey

MADE TO ORDER



Shown in
Waterfall Back

Made
in the UK





Kensey Size Guide:

* Also available in Petite - measurements Seat H 16" x Seat W 18" x Seat D 16" - on Made to order only

	Small	Medium	Large
Seat height	45.7cm / 18"	50.8cm / 20"	53.3cm / 21"
Seat width	45.7cm / 18"	50.8cm / 20"	55.9cm / 22"
Seat depth	45.7cm / 18"	50.8cm / 20"	53.3cm / 21"
Arm height	20.3cm / 8"	19.1cm / 7½"	19.1cm / 7½"
Back height	66cm / 26"	68.6cm / 27"	71.1cm / 28"
Overall width	76.2cm / 30"	81.3cm / 32"	86.4cm / 34"
Overall height	108cm / 42½"	110.5cm / 43½"	113cm / 44½"
Weight limit	20 stone / 127kgs	25 stone / 158.8kgs	25 stone / 158.8kgs

The Kensey is our new UK Premium Tilt-in-Space (TIS) Single-motor Cosi Chair providing a smooth rise and recline movement. It's Waterfall Back design (shown on the front) delivers support to the user's back, shoulders and neck areas throughout travel. Choose from any other back design below:



Main Features:

- Premium Seminar TIS mechanism with Variable Angle Lift
- Arm covers and Antimacassar included Free of Charge as standard
- Choice of any back design or style of fabric on a Made to order basis **only**
- 5-button power on handset
- New pleated wing design
- Heat & massage / drop in seat cushions for pressure care issues - both available on request*
- 5 Year Guarantee on frames, mechanisms and all electrics

Local Cosi Chair Retailer:

*Available on Made to order **only** - please ask for more details.
Please note colours may vary slightly from the images shown.