

SANGO ADVANCED EME

Version: 1.4
 Valid from: 16/04/2018
 Valid until: New version or latest 31-12-2018
 Price list no: 1000010



ELECTRIC MOBILITY

Frontwheel drive: £ 4705,- (excl. VAT)
Rearwheel drive: £ 4705,- (excl. VAT)
Midwheel drive: £ 5585,- (excl. VAT)

Includes the following options

- 6 km/h
- Castor wheels 9", air
- Drive Wheels 14", air
- Suspension on drive wheels
- Frame colours: Red, Blue or Gold
- 50 Ah batteries
- 8A charger
- VR2 70A controller
- Drive only joystick module
- Swing-away joystick bracket, right hand
- Seat height set at 42 cm (min. 39 cm)
- Standard contour seat 46 x 48 cm (w x d)
- Contour standard backrest 46 x 50 cm (w x h)
- Swing away footrests with calf strap
- Flip-up armrests
- PU armpads
- Lapbelt
- Anti-tippers



Type: Order Quotation

Name dealer: _____

Reference no.: _____

Delivery address: _____

Remarks: _____

This picture may deviate from the standard configuration

Specifications

Maximum user weight 160 kg
 Transport brackets ISO7176-19
 CE-marked

Prices in GBP excluding VAT.
 Subject to change without prior notice.

Standard = STD

DRIVE BASE			£
6000001	<input type="checkbox"/>	Front Wheel Drive (FWD)	4705
6000002	<input type="checkbox"/>	Rear Wheel Drive (RWD)	4705
6000003	<input type="checkbox"/>	Mid Wheel Drive (MWD)	5585
Maximum speed			
6000004	<input type="checkbox"/>	6 km/h motor	STD
6000005	<input type="checkbox"/>	6 km/h HT motor	250
6000006	<input type="checkbox"/>	10 km/h motor	(FWD R-net Gyro module mandatory) 275
6000007	<input type="checkbox"/>	10 km/h HT motor	(FWD R-net Gyro module mandatory) 390
6000008	<input type="checkbox"/>	12,5 km/h HT motor	(Only with R-net, 120A PM recommended, max user weight 130 kg)(FWD R-net Gyro module mandatory) 695
Frame colour			
6000010	<input type="checkbox"/>	Red metallic	
6000011	<input type="checkbox"/>	Matt blue metallic	
6000012	<input type="checkbox"/>	Matt Gold metallic	
6000013	<input type="checkbox"/>	Different colour	(please specify RAL colour) 375
Batteries & charger			
6000014	<input type="checkbox"/>	Batteries 50Ah	STD
6000015	<input type="checkbox"/>	Batteries 60Ah	85
6000016	<input type="checkbox"/>	Batteries 78Ah	235
6000017	<input type="checkbox"/>	Charger 24V - 8A	STD
6000018	<input type="checkbox"/>	Charger 24V - 12A	195

WHEELS			£
Wheel size			
6000019	<input type="checkbox"/>	12" Drive wheels with 8" castor wheels	FWD/RWD
6000020	<input type="checkbox"/>	14" Drive wheels with 9" castor wheels	FWD/RWD
6000022	<input type="checkbox"/>	12" Drive wheels with 7" castor wheels	MWD
6000023	<input type="checkbox"/>	14" Drive wheels with 8" castor wheels	MWD
Wheel type			
(air = puncture protected by special layer)			
6000024	<input type="checkbox"/>	Castor wheels air	FWD/RWD
6000025	<input type="checkbox"/>	Castor wheels air	MWD (7" always puncture proof)
6000026	<input type="checkbox"/>	Castor wheels solid	FWD/RWD
6000027	<input type="checkbox"/>	Castor and support wheels solid	MWD
6000028	<input type="checkbox"/>	Drive wheels air	ALL
6000029	<input type="checkbox"/>	Drive wheels solid	ALL
Tyres			
6000030	<input type="checkbox"/>	Black	STD
6000031	<input type="checkbox"/>	Grey	
Mudguards			
6000032	<input type="checkbox"/>	Mudguards castor wheels 9"	41
6000033	<input type="checkbox"/>	Mudguards drive wheels 14"	45

POWERED OPTIONS & LIGHTING			£
Seat adjustments			
6000034	<input type="checkbox"/>	Powered backrest recline 30°	(1 actuator) 425
6000035	<input type="checkbox"/>	Powered seat tilt in space 45°	(1 actuator) 840
6000036	<input type="checkbox"/>	Powered seat lift 30 cm and seat tilt in space 45°	(2 actuators) 1695
Legrest adjustments			
(for manual legrests and legrest accessories, see Manual Legrests)			
6000037	<input type="checkbox"/>	Powered articulating legrests	(6°-70°, angle adjustable footplates, 2 actuators) 831
6000038	<input type="checkbox"/>	Powered centre mount legrest, angle adjustable	(Angle adj restricted for RWD)(length comp., angle adj footplates, 1 actuator) 775
6000039	<input type="checkbox"/>	Powered centre mount legrest, angle and length adj.	(Angle adj restricted for RWD)(length comp., angle adj footplates, 2 actuators) 1075
Lighting			
6000040	<input type="checkbox"/>	LED lights & indicators	Required by UK law if 10km/h or 12,5km/h motors are selected 390

CONTROLS			£
Mounting of the Joystick module			
6000043	<input type="checkbox"/>	Swing-away bracket, right	STD
6000044	<input type="checkbox"/>	Swing-away bracket, left	
Joystick knob options			
6000074	<input type="checkbox"/>	Mushroom knob	50
6000075	<input type="checkbox"/>	Ball mini, Ø 32mm	50
6000076	<input type="checkbox"/>	Softball big, Ø70 mm	50
6000077	<input type="checkbox"/>	Stick joystick	50
6000078	<input type="checkbox"/>	T-grip	59
6000079	<input type="checkbox"/>	U-Grip vertical	75
6000080	<input type="checkbox"/>	U-Grip horizontal	75
VR2 remotes & power modules			
(max 2 actuators)			
Remotes			
6000045	<input type="checkbox"/>	VR2 Drive only	STD
6000046	<input type="checkbox"/>	VR2 Lights	35
6000047	<input type="checkbox"/>	VR2 Seat & lights	60
Power modules			
6000048	<input type="checkbox"/>	70A Power module	STD
6000177	<input type="checkbox"/>	70A Power module, 2 actuators	65
6000049	<input type="checkbox"/>	90A Power module, 2 actuators	130
Attendant control			
6000050	<input type="checkbox"/>	VR2 + attendant control, right	(dual) 450
6000051	<input type="checkbox"/>	VR2 + attendant control, left	(dual) 450

R-net remotes, power modules and options			
Remotes (all underneath LED and LCD joystick modules can operate 6 actuators and lights)			
6000052	<input type="checkbox"/>	R-net LED	(If no power module selected, 90A will be chosen) 0
6000053	<input type="checkbox"/>	R-net LCD-1	(CJSM1) (2 jack sockets) 140
6000055	<input type="checkbox"/>	R-net LCD-2, Infra-red	(CJSM2)(2 jack sockets) 395
6000057	<input type="checkbox"/>	R-net LCD-2, Infra-red and Bluetooth	(CJSM2)(2 jack sockets) 635
6000059	<input type="checkbox"/>	R-net LCD-1, for tray control (please choose also tray)	(CJSM1) 550
Power modules			
6000060	<input type="checkbox"/>	90A Power module	465
6000061	<input type="checkbox"/>	120A Power module	(recommended with 12,5 km/h) 580
6000062	<input type="checkbox"/>	DAC box, direct actuator control	(only with R-net LCD joystick modules) 350
6000063	<input type="checkbox"/>	Gyro module	(mandatory on FWD 10 and 12,5 km/h) 295
6000064	<input type="checkbox"/>	R-net Bluetooth module	(for Android and Windows) 625
6000065	<input type="checkbox"/>	R-net i-Device Bluetooth module	(for Apple IOS) 625
6000066	<input type="checkbox"/>	R-net 2nd Bluetooth module	(for Android and Windows) 625
Attendant control			
6000067	<input type="checkbox"/>	R-net + attendant control, right	(dual) 450
6000068	<input type="checkbox"/>	R-net + attendant control, left	(dual) 450
Omni module (Special controls, see end of order form)			
6000070	<input type="checkbox"/>	Mounted on fixed bracket, right	825
6000071	<input type="checkbox"/>	Mounted on fixedbracket, left	825
6000072	<input type="checkbox"/>	Mounted on Swing-away bracket, right	875
6000073	<input type="checkbox"/>	Mounted on Swing-away bracket, left	875
Push button switches			
6000081	<input type="checkbox"/>	Small, red	90
6000082	<input type="checkbox"/>	Small, green	90
6000083	<input type="checkbox"/>	Big, red	90
6000084	<input type="checkbox"/>	Big ,green	90

SEATS			£
Seat height (floor to top seat frame in cm with 14" drive wheels and 0° tilt)			
6000085	<input type="checkbox"/>	Position -1*: SH 39 cm, SD <52 cm	
6000086	<input type="checkbox"/>	Position -1*: SH 39 cm, SD >52 cm	
6000087	<input type="checkbox"/>	Position 0: SH 42 cm	STD
6000088	<input type="checkbox"/>	Position 1: SH 44 cm	
6000089	<input type="checkbox"/>	Position 2: SH 46 cm	
6000090	<input type="checkbox"/>	Position 3: SH 48 cm	
Seat cushion size			
6000184	<input type="checkbox"/>	Width 42 cm	
6000185	<input type="checkbox"/>	Width 44 cm	
6000186	<input type="checkbox"/>	Width 46 cm	STD
6000187	<input type="checkbox"/>	Width 48 cm	
6000188	<input type="checkbox"/>	Width 50 cm	
6000189	<input type="checkbox"/>	Width 52 cm	
6000190	<input type="checkbox"/>	Width 54 cm	99
6000191	<input type="checkbox"/>	Depth 48 cm	STD
6000192	<input type="checkbox"/>	Depth 55 cm	19
6000091	<input type="checkbox"/>	Different seat cushion size (specify desired size)	215
Seat cushion type			
6000092	<input type="checkbox"/>	Contour	finish: vinyl STD
6000094	<input type="checkbox"/>	Contour, partially strapped	air mesh 95
6000096	<input type="checkbox"/>	Anatomic Formed	air mesh (partially strapped) 170
6000098	<input type="checkbox"/>	Anatomic Formed Visco	air mesh (partially memory foam and webbed) 195

	Seat boards	finish:	
6000100	<input type="checkbox"/> Seat board	vinyl	0
6000102	<input type="checkbox"/> Seat board, partially strapped	vinyl	95
6000104	<input type="checkbox"/> Threshold, front	vinyl	(upholstered, H=5cm) 35
6000105	<input type="checkbox"/> Threshold, rear	vinyl	(upholstered, H=5cm) 35
6000106	<input type="checkbox"/> Acrylic threshold, rear	acryl	(depth adjustable, H=5cm) 39
	Seat cushion finish		
6000107	<input type="checkbox"/> Air mesh instead of vinyl		75
6000108	<input type="checkbox"/> Stretchrubber instead of air		40
	Seat accessories		
6000110	<input type="checkbox"/> Thigh support small, set		(22 x 11cm LxH) 275
6000111	<input type="checkbox"/> Thigh support large, set		(36 x 13cm LxH) 310
6000112	<input type="checkbox"/> Multi thigh support small, set		(22 x 11cm LxH) 375
6000113	<input type="checkbox"/> Multi thigh support large, set		(36 x 13cm LxH) 410
6000232	<input type="checkbox"/> Acrylic shield for (multi adjustable) thigh support instead of upholstered		0

BACKRESTS			£
	Backrest cushion size		
6000193	<input type="checkbox"/> Width 42 cm		
6000194	<input type="checkbox"/> Width 44 cm		
6000195	<input type="checkbox"/> Width 46 cm		STD
6000196	<input type="checkbox"/> Width 48 cm		
6000197	<input type="checkbox"/> Width 50 cm		
6000198	<input type="checkbox"/> Width 52 cm		
6000199	<input type="checkbox"/> Width 54 cm		99
6000200	<input type="checkbox"/> Height 50 cm		STD
6000201	<input type="checkbox"/> Height 57 cm		19
6000115	<input type="checkbox"/> Different backrest cushion size (specify desired size)		225
	Backrest cushion type	finish	
6000116	<input type="checkbox"/> Contour	vinyl	STD
6000117	<input type="checkbox"/> Deep contour	air mesh	(width adj. built-in laterals) 110
6000118	<input type="checkbox"/> Deep contour vertical	air mesh	(vertical elastic straps and width adj. built-in laterals) 175
6000119	<input type="checkbox"/> Deep contour horizontal	air mesh	(horizontal tension adj. straps and width adj. built-in laterals) 195
6000120	<input type="checkbox"/> Deep contour Flex	air mesh	(2 wooden strips, horiz.tension adj. straps and width adj. built-in laterals) 210
6000121	<input type="checkbox"/> Extra Deep contour	air mesh	(horizontal tension adj. straps and width adj. built-in laterals) 255
6000122	<input type="checkbox"/> Extra Deep contour Flex ,	air mesh	(2 wooden strips, horiz.tension adj. straps and width adj. built-in laterals) 270
	Backrest cushion finish		
6000123	<input type="checkbox"/> Air mesh instead of vinyl		75
6000124	<input type="checkbox"/> Stretchrubber instead of air		40
6000125	<input type="checkbox"/> Different backrest cushion finish (specify desired finish)		225
	Backrest accessories		
6000126	<input type="checkbox"/> Thoracic support, adjustable, right		275
6000127	<input type="checkbox"/> Thoracic support, adjustable, left		275
6000128	<input type="checkbox"/> Manual adjustable recline (gas strut)		210
6000129	<input type="checkbox"/> Folding backrest		125

MANUAL LEGRESTS			£
		(for powered legrests see Powered options and lighting)	
	Centre mount legrest		
6000130	<input type="checkbox"/> Centre mount legrest, angle and height adjustable		0
	Centre mount legrest accessories		
6000131	<input type="checkbox"/> Calf plates PU		145
	Swing-away legrests		
6000132	<input type="checkbox"/> Standard footrests	(angle adjustable foot plates and calf strap)	STD
6000133	<input type="checkbox"/> Standard footrests with PU calf plates	(depth adj. calf plates and angle adjustable foot plates)	145
6000134	<input type="checkbox"/> Comfort legrests	(6°-70° gas spring operated , depth adjustable PU calf plates and angle adjustable foot plates)	275
	Swing-away legrest accessories		
6000135	<input type="checkbox"/> Calf plates upholstered instead of PU		48
6000136	<input type="checkbox"/> Kneepad, right		46
6000137	<input type="checkbox"/> Kneepad, left		46
6000138	<input type="checkbox"/> Canvas heel support, right		44
6000139	<input type="checkbox"/> Canvas heel support, left		44
6000140	<input type="checkbox"/> Toe strap, right		36
6000141	<input type="checkbox"/> Toe strap, left		36
6000142	<input type="checkbox"/> Amputation support, right		127
6000143	<input type="checkbox"/> Amputation support, left		127

ARMRESTS			£
	Armrests		

6000145	<input type="checkbox"/>	Flip-up	(angle- and height adjustable)	STD
Armpads				
6000146	<input type="checkbox"/>	PU 34 cm, right		STD
6000147	<input type="checkbox"/>	PU 34 cm, left		STD
6000148	<input type="checkbox"/>	Upholstered 25 cm, right		55
6000149	<input type="checkbox"/>	Upholstered 25 cm, left		55
6000150	<input type="checkbox"/>	Upholstered 40 cm, right		60
6000151	<input type="checkbox"/>	Upholstered 40 cm, left		60

HEADRESTS				£
Headrest metal bracket				
6000154	<input type="checkbox"/>	Adjustable in height and depth		90
6000155	<input type="checkbox"/>	Adjustable in height, depth and width		116
6000156	<input type="checkbox"/>	Quick release levers		29
Headrest cushion				
6000157	<input type="checkbox"/>	PU small		90
6000158	<input type="checkbox"/>	Padded butterfly shape		99
6000159	<input type="checkbox"/>	Padded flexible small		115
6000160	<input type="checkbox"/>	Padded flexible large		130

TRAYS				£
Acrylic tray				
6000161	<input type="checkbox"/>	Side joystick recess	(mounted at opposite side to the joystick module)	255
6000162	<input type="checkbox"/>	Hemi tray	(half tray)(mounted at opposite side to the joystick module)	245
6000163	<input type="checkbox"/>	Joystick recess in centre, tube mounted right		275
6000164	<input type="checkbox"/>	Joystick recess in centre, tube mounted left		275
Tray accessories				
6000165	<input type="checkbox"/>	Partially upholstered		85
6000166	<input type="checkbox"/>	Raised edge acrylic, front		42
6000167	<input type="checkbox"/>	Raised edge acrylic, right		42
6000168	<input type="checkbox"/>	Raised edge acrylic, left		42

ACCESSORIES				£
6000169	<input type="checkbox"/>	Lapbelt standard, 1.90 m		STD
6000170	<input type="checkbox"/>	Lapbelt long, 2.90 m		35
6000298	<input type="checkbox"/>	Mirror on right hand side	Is required by UK law if 10km/h or 12,5km/h motors are selected	41
6000171	<input type="checkbox"/>	Mirror on joystick side		41
6000172	<input type="checkbox"/>	Mirror not on joystick side		43
6000173	<input type="checkbox"/>	Cane/crutch holder for two walking canes		59
6000174	<input type="checkbox"/>	Kerb climber	(only with RWD, not with centre mount legrest)	210
6000175	<input type="checkbox"/>	Dahl car docking preparation		95

Special controls

Different delivery time

Additional on chair, as spare part see spare part list

Prices in GBP excluding VAT.

Subject to change without prior notice.

Standard = STD

SPECIAL CONTROLS		£
OMNI Special input devices		
		(operating force)
8000003	<input type="checkbox"/> Micro joystick	(8,5 gram) 2357
8000004	<input type="checkbox"/> Multi joystick	(50 gram) 2130
8000005	<input type="checkbox"/> All-round Light joystick	(120 gram) 1694
8000006	<input type="checkbox"/> All-round Light joystick , chin	(incl. 2 satellite switches en foam ball)(120 gram) 2058
8000007	<input type="checkbox"/> All-round joystick	(250 gram) 1065
8000008	<input type="checkbox"/> All-round joystick , chin	(incl. 2 satellite switches en foam ball)(250 gram) 1427
8000031	<input type="checkbox"/> All-round heavy duty joystick	(650 gram) 2308
R-net (direct plug-in) Special input devices		
		(operating force)
8000032	<input type="checkbox"/> Micro joystick	(8,5 gram) 2475
8000033	<input type="checkbox"/> Multi joystick	(50 gram) 2248
8000034	<input type="checkbox"/> All-round Light joystick	(120 gram) 1810
8000035	<input type="checkbox"/> All-round Light joystick , chin	(incl. 2 satellite switches en foam ball)(120 gram) 2176
8000036	<input type="checkbox"/> All-round joystick	(250 gram) 1186
8000037	<input type="checkbox"/> All-round joystick , chin	(incl. 2 satellite switches en foam ball)(250 gram) 1549
8000038	<input type="checkbox"/> All-round heavy duty joystick	(650 gram) 2427
Mounting of special input device		
8000009	<input type="checkbox"/> On armrest, right	84
8000010	<input type="checkbox"/> On armrest, left	84
8000011	<input type="checkbox"/> In /on tray	735
		specify dimensions and required mounting position in a detailed drawing
Swing-away arms for chin control:		
8000012	<input type="checkbox"/> Manual horizontal, right	1337
8000030	<input type="checkbox"/> Manual horizontal, left	1337
8000013	<input type="checkbox"/> Powered horizontal, right	2487
8000014	<input type="checkbox"/> Powered horizontal, left	2487
8000015	<input type="checkbox"/> Powered vertical , right	2487
8000016	<input type="checkbox"/> Powered vertical, left	2487
8000039	<input type="checkbox"/> Switch Ø 17mm to operate swingarm, right	(mounted on headrest bracket) 427
8000040	<input type="checkbox"/> Switch Ø 17mm to operate swingarm, left	(mounted on headrest bracket) 427
Switches		
	<input type="checkbox"/> Please specify required type, location and way of mounting of the switches	POR
REMARKS:		

