

The

Ellen

Tilt-in-Space Single-motor

Made
in the UK



Ellen

Waterfall Back

Choose our furniture with confidence -
All riser recliner chairs are tested by FIRA -
The leading authority in furniture safety by
design.



Riser Recliners available from STOCK in the following fabrics



Seasons
Mink



Seasons
Oatmeal



Seasons
Plum



Seasons
Cocoa



2 or 3-Seater static sofa available
on Made to order

Colours may vary slightly from images shown which are for illustrative purposes only.

Main Features:

The Ellen Cosi Chair is a comfortable Tilt-in-Space Single-motor chair with pocket sprung chaise seating and a deep fibre filled Waterfall backrest. The design retains the same angle between the seat and backrest to maximise comfort for individuals with delicate skin conditions or where legs need elevating above the hips.

- 3-button power on handset
- Waterfall Back design
- Pocket sprung seat base
- Scroll arms with wings
- Arm covers and Antimacassar included
- Other fabrics are available on Made to order
- 2 or 3-Seater static sofa available on Made to order
- 5 Year Guarantee on frames, mechanisms & all electrics
- Quality Craftsmanship, **BRITISH MADE**

Ellen Size Guide:

	One size
Seat height	48.3cm / 19"
Seat width	50.8cm / 20"
Seat depth	50.8cm / 20"
Arm height	20.3cm / 8"
Back height	71.1cm / 28"
Overall width	82.5cm / 32½"
Overall height	111.8cm / 44"
Weight limit	20 stone / 127kg

Local Cosi Chair Retailer: