

The

Medina

Independent Dual-motor

Made
in the UK



CosiChair

Medina

Waterfall Back

Choose our furniture with confidence -
All riser recliner chairs are tested by FIRA -
The leading authority in furniture safety by
design.



Riser Recliners available from STOCK in the following fabrics



Colours may vary slightly from images shown which are for illustrative purposes only.

The Medina Cosi Chair is a superbly comfortable independent Dual-motor chair with pocket sprung chaise seating and a soft fibre filled backrest. The independent control of backrest and footrest functions operated through the simple handset enables any individual to find their desired seating or reclined position for maximum comfort.

Main Features:

- 4-button power on handset
- Waterfall, Lateral or Button Back design
- Pocket sprung seat base
- Scroll arms with wings
- Arm covers and Antimacassar
- Other fabrics are available on Made to order
- 2 or 3-Seater sofa available on Made to order
- 5 Year Guarantee on frames, mechanisms and all electrics
- Quality Craftsmanship, **BRITISH MADE** 

Local Cosi Chair Retailer:

Medina Waterfall Back Size Guide:

	Small	Medium	Large
Seat height	46.7cm / 18½"	49.5cm / 19½"	53.3cm / 21"
Seat width	46.7cm / 18½"	50.8cm / 20"	55.9cm / 22"
Seat depth	45.7cm / 18"	50.8cm / 20"	53.3cm / 21"
Arm height	19.1cm / 7½"	19.1cm / 7½"	19.1cm / 7½"
Back height	66cm / 26"	68.6cm / 27"	71.1cm / 28"
Overall width	81.3cm / 32"	83.8cm / 33"	86.4cm / 34"
Overall height	108cm / 42½"	110.5cm / 43½"	113cm / 44½"
Weight limit	25 stone / 158.8kg	25 stone / 158.8kg	25 stone / 158.8kg

