

Rivco

Mid-wheel powerchairs



Your new route to independence

RASCAL 



Rivco Specification:

Type Class	B / 2
Overall length	110cm / 43.3"
Overall width	70cm / 27.5"
Floor to seat height	54-56.5cm (min-max)
Seat width	18.5", 20.5", 22"
Seat depth	18.9", 19.5", 19.5"
Maximum carrying capacity	204kgs / 32 stone
Battery voltage & capacity	12V / 55Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	5cm / 2"
Maximum kerb height	5cm / 2"
Range*	up to 30.2km / 18.8miles
The maximum safe slope	6°
Turn-around width	1.60 / 5.2'
Colour	Blue Moon / Sunset Orange

*Range on full charge and flat ground
See full Tech Page on Website & Owner's Manual for details

Rascal Rivco Profile:

The Rascal Rivco is a heavy-duty 4 mph Mid-Wheel Drive outdoor/indoor powerchair, with durable full suspension and comfortable cushioned Captain's seat mounted on a larger base. The Rivco delivers great performance with its 55Ah batteries contributing to its excellent range, outstanding manoeuvrability and a smooth ride.

Ryley, Ryley SL and Rivco Seat Lift models also available.

Local Rascal Mobility Retailer:

Main Features:

- Flip-up fully adjustable armrests
- Full suspension
- Height adjustable Swing-away joystick
- Convenient joystick charging port
- Adjustable Captain's seat
- Flip-up footplate

OPTIONAL UPGRADES

Optional Upgrades:

- Swing-away legrests
- Elevating legrests
- Attendant Control

OPTIONAL ACCESSORY

Optional Extras:

- Lighting Kit
- Armrest bag
- Rear Saddle bag
- Cane & Crutch holder bag
- Oxygen bottle bag



Specifications subject to change without prior notice
Images are for illustrative purposes only
E. and O.E- v 1 Final-January 2023