

Veo Sport **LiFe**

Lithium powered scooters

Great Features



LiFetime Guarantee

on Lithium Ferro-Phosphate Battery*



The Lithium Ferro-Phosphate Battery is **'built to last'** & the **SAFEST** Lithium Battery in the Mobility sector



Veo Sport **LiFe**

LiFetime Guarantee
on Lithium Ferro-Phosphate Battery†

†T&Cs apply



Veo Sport **LiFe** Specification:

Type Class	B / 2
Overall length	108cm / 42.5"
Overall width	60cm / 23.6"
Floor to seat height	47-49cm (min-max)
Seat width	45.7cm / 18"
Seat depth	40.6cm / 16"
Maximum carrying capacity	130kgs / 20.5 stone
Battery voltage and capacity	25.6V / 6.2Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	6cm / 2.4"
Maximum kerb height	5cm / 2"
Range*	up to 23.2km / 14.4miles
The maximum safe slope	6°
Turn-around width	1.40m / 4'6"
Colour	Neon Blue / Mulberry
Single Heaviest component	20.2kgs / 44.5lbs

*Range on full charge and flat ground
See full Tech Page on Website & Owner's Manual for details

Rascal Veo Sport **LiFe** Profile:

Offering peace of mind with a LiFetime Guarantee on the battery, it's all about the State-of-The-Art Lithium Ferro-Phosphate Battery weighing just 4.7kg on this striking Rascal Veo Sport LiFe. Loaded with impressive features & modern, sporty design with class leading all-round suspension. With an Airline Compliant 12.4Ah battery, offering a range of up to 14.4 miles, plus padded delta tiller boasting a convenient charging point and luxury seat, detailed with unique white accent stitching, the Veo Sport LiFe is a real pleasure to ride!

The Veo Sport LiFe's easy and robust take apart mechanism allows you to store or drive away with it in the boot of your car, so it is very practical yet with its distinctive looks, exclusive all-round suspension plus LED front light, this is the travel scooter with the edge.

Air travel & CAA Lithium Battery Legislation



Airline Compliant 12.4Ah battery meets I.A.T.A. UN38.3 Passenger Aircraft Transportation requirements for overseas travel & can be taken on a plane.

(Subject to Airline Regulation & Restrictions)


† LiFetime Battery Guarantee T&Cs apply



Veo Sport SR Model with standard 22Ah or lighter 14Ah Battery options also available

Main Features:

- British Made / Medical /Military grade **Lithium Ferro-Phosphate Battery**
- Car transportable
- All-round suspension
- Stylish sports seat
- **Tiller charging point**
- **LED front light**
- **Front bumper**
- Dual compound puncture proof tyres
- Comfortable delta tiller



12.4Ah weighs just 4.7kg:

Relax & enjoy the ride in the knowledge that the **Lithium Ferro-Phosphate** (LiFePO4) Battery pack is the safest, **greenest** & is warranted for a design life of over 1500 cycles providing an increased efficiency & performance to stand up to everyday life.

The LED front light & bumper together with class leading all-round suspension unites design, functionality & safety.



Specifications subject to change without prior notice

Images are for illustrative purposes only
E. and O.E- v 2 Aug 2023

electricmobility.co.uk