

Vie Life

Lithium powered scooters

LiFetime Guarantee
on Lithium Ferro-Phosphate Battery*



RASCAL
Life **BRITISH**
MADE



The Lithium Ferro-Phosphate Battery is **'built to last'** & the **SAFEST** Lithium Battery in the Mobility sector

RASCAL



Vie LiFe 10Ah Specification:

Type Class	B / 2
Overall length	108.3cm / 42.6"
Overall width	47cm / 18.6"
Floor to seat height	39-43cm (min-max)
Seat width	41cm / 16.1"
Seat depth	37cm / 14.6"
Maximum carrying capacity	121kgs / 19 stone
Battery voltage & capacity	25.6V / 10Ah
Maximum speed	6.4km/h / 4mph
Ground clearance	4.5cm / 1.8"
Maximum kerb height	5cm / 2"
Range*	up to 19.2km / 11.9 miles
The maximum safe slope	6°
Turn-around width	1.25m / 4'1"
Colour	Blue / Red
Single Heaviest component	14.8kgs / 32.5lbs

*Range on full charge and flat ground

See full Tech Page on Website & Owner's Manual for details

Rascal Vie LiFe Profile:

Here comes affordability for everyone. Economically priced, with a range upto 11.9 miles makes it suitable for your daily run about. Offering peace of mind with a LiFetime Guarantee on the battery, it's all about the State-of-The-Art Lithium Ferro-Phosphate Battery weighing just 3.7kg. This 4mph scooter is comfortable to use, yet rugged & nimble with a tight turning circle. Transportable with the front section weighing just 14.8kg, giving you all round versatility & daily enjoyment when you are outdoors with friends & family. The Rascal Vie LiFe is a very practical scooter with an easy & robust take apart mechanism which allows you to store or drive away with it in the boot of your car, perfect for everyday. Experience this class-leading boot scooter for yourself.

Air travel & CAA Lithium Battery Legislation




Airline Compliant 10Ah battery meets I.A.T.A. UN38.3 Passenger Aircraft Transportation requirements for overseas travel & **can be taken on a plane.**

(Subject to Airline Regulation & Restrictions)

† LiFetime Battery Guarantee T&Cs apply

Main Features:

- British Made / Medical / Military grade **Lithium Ferro-Phosphate Battery**
- **A quality product, economically priced**
- Manoeuvrable with tight turning circle
- Car transportable
- Delta tiller & convenient charging point
- Comfort swivel seat with white stitching
- Solid foam filled puncture proof tyres
- Spacious shopping basket
- Rear safety reflectors



10 Ah weighs just 3.7kg:

Relax & enjoy the ride in the knowledge that the **Lithium Ferro-Phosphate (LiFePO4)** Battery pack is the safest, **greenest** with an increased efficiency & performance to stand up to everyday life.



Specifications subject to change without prior notice
Images are for illustrative purposes only
E. and O. E- v1 August 2024