

# Vista **Life**

Lithium powered scooters

**LiFetime Guarantee**

on Lithium Ferro-Phosphate Battery\*



**RASCAL**  
**Life** **BRITISH**  
**MADE**



The Lithium Ferro-Phosphate Battery is **'built to last'** & the **SAFEST** Lithium Battery in the Mobility sector

**RASCAL**



## Vista LiFe 36Ah Specification

Type Class	B / 2
Overall length	118cm / 46.5
Overall width	58.5cm / 23"
Floor to seat height	57cm - 62cm
Seat width	44.5cm / 17.5"
Seat depth	43.1cm / 17"
Maximum carrying capacity	150kgs / 23.6 stone
Battery voltage and capacity	12.8V / 36Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	7.5cm / 3"
Maximum kerb height	6cm / 2.3"
Range*	up to 53.9km / 33.5 miles
The maximum safe slope	6°
Turn-around width	1.7m / 5'8"
Colour	Oxford Blue / Sunset Orange

Single Heaviest component (front Chassis) 25kgs / 55.12 lbs

\*Range on full charge and flat ground

See full Tech Page on Website & Owner's Manual for details

## Main Features:



- British Made / Medical / Military grade 36Ah **Lithium Ferro-Phosphate Battery**
- All-round suspension
- Easy take apart action - car transportable
- Comfort seat with folding backrest, adjustable arms & integrated swivel function
- Delta tiller with useful charging port
- Class leading weight capacity
- Off-Board Battery Charging capability

## Alternative DX option:

Featuring:

- Captain's seat with folding backrest, headrest, adjustable arms with swivel & slide functions
- Seat width: 45.7cm / 18"
- Seat depth: 43.1cm / 17"



## Rascal Vista **LiFe** Profile:

Introducing the Rascal Vista LiFe with a LiFetime Battery Guarantee. This mid-size pavement scooter benefits from the quick release feature, it dismantles into a few easy-to-lift sections for transporting in your car. The compact & slim-line style offers plenty of legroom yet it is still extremely robust with generous turning manoeuvrability, making it ideal for everyday use. Available in two eye catching colours with appealing features including easy access tiller charging & delta handle bar plus the all-round suspension offers added comfort. **When taking apart, the split battery pack is also a real advantage.** It comes with front & rear LED lights maximising visibility, giving you both confidence & support wherever you go. Available with solid puncture proof or pneumatic tyres - the choice is yours!



36Ah split battery pack, weighing just 5.35kg each:

Relax & enjoy the ride in the knowledge that the **Lithium Ferro-Phosphate** (LiFePO<sub>4</sub>) Battery pack is the safest, **greenest** with an increased efficiency & performance to stand up to everyday life. (Total weight of battery 10.7kg)



† LiFetime Battery Guarantee T&Cs apply



Vista model with 35Ah & Vista DX model with 50Ah battery options also available

Specifications subject to change without prior notice  
Images are for illustrative purposes only  
E. and O.E- v2 September 2024